


System knee joint set with lock Artikon.genu-bi-lock

High-quality system joints for the construction of bilateral orthoses: The new monocentric Artikon.genu-bi-lock knee joints have a manual locking mechanism that can be released by the wearer either with the attached locking lever or an optional release device.

Characteristics

- system knee joint with locking mechanism for movement guidance
- with 0° extension locked
- unlocking via lever or release cable 2J35
- optimum positioning of the rails by repositioning the joints at the pivot point
- with slide bearing





Streifeneder
ortho.production GmbH
Moosfeldstrasse 10
82275 Emmering
Germany



T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W151/E

