

Preview: Technology can be so beautiful – new design for KINEGEN.stream


The polycentric prosthetic knee joint KINEGEN.stream presents itself in new product design. The knee joint's main colour is a classy titan grey, which provides a modern look in combination with white coloured areas. Hydraulic components in glowing red enhance the lively look. The dynamic-organic shape with soft edges and smooth feel underline the natural, positive impression.

Technical data

- max. body weight: 150 kg
- activity level (walky): 3 to 4
- effective height: 200 mm

New design available from 4th quarter of 2016!





Streifeneder
ortho.production GmbH
Moosfeldstraße 10
82275 Emmering
Germany



T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W130/E