



Design without sharp edges: Pre-preg cord 150C1

The pre-preg cord allows an even edge design of orthoses and prevents sharp rimlines. This optimises wearing comfort for the user and prevents flexible inlays or clothing from being worn through.

Characteristics

- available in thicknesses of 1 mm or 3 mm
- BBD as with regular pre-preg
- cord is incorporated into the bidirectional pre-preg
- cord is protected by peel-off foil
- it can be cut to exact lengths





Streifeneder
ortho.production GmbH
Moosfeldstrasse 10
82275 Emmering
Germany



T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W155/E

