


## Deep-drawing in one step: Streifylit Duo 111P500

The new Streifylit Duo convinces by its unique material combination of plastic and PO-film (foam cover). This material enables to manufacture quickly and efficiently individual, padded night splints in just one deep-drawing process.

### Characteristics

- optimal processing temperature at approx. 160°C / 320°F
- for deep-drawing, colourable, light weight, low shrinkage
- break- and scratch resistant, easy to cut, grind and polish
- bio compatible, disinfectable
- PU = 1 sheet (1000 x 1200 x 4,5 mm)





Streifeneder  
ortho.production GmbH  
Moosfeldstraße 10  
82275 Emmering  
Germany



T +49 8141 6106-0  
F +49 8141 6106-70  
export@streifeneder.de  
www.streifeneder.com/op

300W159/E