

PediSupport.dorsal Secure support for the foot and heel



© Africa Studio / fotolia

Start-up pains in the morning or after considerable periods of inactivity, including piercing sensations in the heel area can be symptoms indicating a condition known as plantar fasciitis. These complaints are often associated with a heel spur. PediSupport.dorsal (Item-No. 191R1) helps in treating this ailment by slightly stretching the plantar fascia and thus ameliorating the pain. The orthosis is worn at night and it keeps the patient's foot at a 90° position. PediSupport.dorsal is hence used even in talipes equinus positions so as to prevent shortening of the Achilles tendon.

Indications

- Plantar fasciitis
- Inflammation of the Achilles tendon
- Talipes equinus position
- Metatarsalgia
- for night resting
- for contracture prophylaxis

Size	Shoe size	Colour	Item-No.
M	36 – 39	black	191R1/M
L	40 – 44	black	191R1/L

Further therapy options for splayfoot contracture: PediSupport.night (190R1/...)

