

Socket material Streifyflex Black



Streifyflex Black

For flexible innersockets, especially designed for the above knee area (ISNY-technology). Also applicable for arm prosthetics.

Streifyflex is a permanent elastic material, suitable for deep drawing, very flexible, disinfected and thermomoldable. After material processing, visible shrinking traces are not recognisable.

Moldable at approximately 130 – 135 °C.

Item-No.	Dimension L x W x T	PU
111P71/9	400 x 400 x 9 mm	sheet
111P71/12	400 x 400 x 12 mm	sheet
111P71/14	400 x 400 x 14 mm	sheet
111P71/18	400 x 400 x 18 mm	sheet

